


# INTEROSS<sup>®</sup> SYRINGE

ANORGANIC CANCELLOUS  
BONE GRAFT GRANULES



## WHAT IS INTEROSS<sup>®</sup> SYRINGE

InterOss<sup>®</sup> syringe is a natural hydroxyapatite bone grafting material for bone regenerative application in dentistry. Made from a proven multi-step purification process, InterOss<sup>®</sup> is an osteoconductive material for bone regeneration. InterOss<sup>®</sup> syringe is designed to deliver InterOss<sup>®</sup> granules more precisely to the intended treatment site without having to use sterile instruments. The syringe is supplied sterile and is dedicated for single use.



InterOss<sup>®</sup> Syringe Type

Size	0.25-1.0 mm		
Reference	IOSGS025	IOSGS050	IOSGS100
Volume	0.25cc	0.50cc	1.00cc

Size	1-2 mm	
Reference	IOLGS050	IOLGS100
Volume	0.50cc	1.00cc



# SIGMAGRAFT


## THE BONE GRAFT SUBSTITUTE SPECIALISTS

### INDICATIONS FOR USE


InterOss<sup>®</sup> syringe is recommended for:

- ✔ Augmentation or reconstructive treatment of the alveolar ridge
- ✔ Filling of periodontal defects
- ✔ Filling of defects after root resection, apicoectomy, and cysectomy
- ✔ Filling of extraction sockets to enhance preservation of the alveolar ridge
- ✔ Elevation of the maxillary sinus floor
- ✔ Filling of periodontal defects in conjunction with products intended for Guided Tissue Regeneration (GTR) & Guided Bone Regeneration (GBR)
- ✔ Filling of peri-implant defects in conjunction with products intended for GBR


### BENEFITS



Micropore (<10µm)



Mesopore (10-100µm)



Macropore (>100µm)

- ✔ Simple to Use, clean handling, time saving, and convenient
- ✔ Biocompatible and osteoconductive making it ideal for bone regeneration uses
- ✔ Hydrophilic properties allow rapid absorption of patient's blood making it easy to apply
- ✔ Interconnected network of macro and micro pores providing an ideal environment for cell attachment and bone growth

**Contact Us**

Address: SigmaGraft, Inc. 575 Sally Place Fullerton, CA 92831, USA  
Phone: (714) 525-0114, Toll Free: (888) 499-0114  
E-Mail: [info@sigmagraft.com](mailto:info@sigmagraft.com)